


CHIROPRACTIC COUCH · ELECTRIC

ALMA PRO

Ref. EL07BPRO23SL.2.A66.GF



Dimensional Drawing



Newly designed chiropractic table with a wide range of functions – for drop, flexion and distraction techniques. Adjustable head sections increase versatility; tension controls and a removable T-bar flexion control allow free movement around the table.

PRODUCT

Manufacturer	COINFY CARE
Reference	EL07BPRO23SL.2.A66.GF
Category	Chiropractic manipulation table
Version	Electric (Hi-Lo)

DIMENSIONS

Overall length	185 – 200 cm
Width	65 / 55 cm
Height adjustment	56 – 90 cm

DRIVE & OPERATION

Height adjustment	electric (Hi-Lo), foot bar
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER

Padding thickness	7.5 cm
Foam density	100 kg/m³

PERFORMANCE & COMPLIANCE

Warranty	2 years on frame, electrical components and upholstery (except in case of misuse)
----------	---

Features & Equipment

- Electric height (Hi-Lo) and traction operation via foot bar and button
- Head section: width-adjustable padded platform, gas-spring headrest ($\pm 20^\circ$), 10 cm height adjustment, lateral movement $\pm 15^\circ$
- Auxiliary straps included
- Four drop sections: cervical, thoracic, lumbar and pelvic
- Motorised adjustable lumbar traction mechanism
- Abdominal separation: thoracic section can be raised for more diverse treatment approaches
- Vertical tension control and horizontal flexion of the lower body via T-bar
- 4 retractable castors with central locking for easy transport

SHIPPING & PACKAGING Packaging size 176 × 72 × 74 cm · Shipping weight approx. 118 kg · Shipping volume 0.94 m³