


CHIROPRACTIC COUCH · ELECTRIC

APPA

Ref. ELT08F.1.A66.WF



Dimensional Drawing



High-quality chiropractic treatment table with centre column. 6 sections and electric height adjustment via foot switch. 4 pedal-adjustable drops: cervical, thoracic, lumbar and pelvic. Durable, antibacterial upholstery. 2-year warranty.

PRODUCT

Manufacturer	COINFYCARE
Reference	ELT08F.1.A66.WF
Category	Chiropractic manipulation table
Version	Electric (Hi-Lo)

DIMENSIONS

Overall length	185 – 198 cm
Width	62 / 55 cm
Height adjustment	60 – 80 cm

DRIVE & OPERATION

Height adjustment	electric (Hi-Lo), foot bar
Motor	TiMOTION, high reliability
Lifting time	23 seconds

UPHOLSTERY & COVER

Padding thickness	7 cm
Foam density	100 kg/m³
Cover colours	blue or dark grey

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	2 years on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric height adjustment via foot switch
- Headrest adjustable in height and opening
- Extendable footrest
- 4 pedal-adjustable drops: cervical, thoracic, lumbar and pelvic
- Elegant design with two unidirectional castors
- Fire-retardant, highly wear-resistant vinyl upholstery
- Paper-roll holder included

SHIPPING & PACKAGING Packaging size 201 × 73 × 82 cm · Shipping weight approx. 152 kg · Shipping volume 1.2 m³