


EXAMINATION & TREATMENT COUCH · ELECTRIC

BATEC

Ref. 75.EL02C.2.A66.WF



Dimensional Drawing



Robust two-section cardiology table for high load capacity and patient comfort. Upholstery with stitched edges, centrally lockable castors, electric backrest and two cardiac cut-outs for optimal access. Standard foot switch – ready to use without assembly.

PRODUCT

Manufacturer	COINFYCARE
Reference	75.EL02C.2.A66.WF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	195 × 70 cm
Height adjustment	58 – 102 cm
Backrest angle	+75°

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
-------------------	---------------------------------

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³

CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016

Features & Equipment

- Electric backrest for smooth, effortless operation
- Two cardiac cut-outs for easier patient access
- Comfortable upholstery with stitched edges for greater durability
- Centrally lockable castors
- Standard foot switch – ready to use without assembly
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 199 × 75 × 61 cm · Shipping weight approx. 65 kg · Shipping volume 0.91 m³