


EXAMINATION & TREATMENT COUCH · ELECTRIC

BOBATH

Ref. EL02W.1.A66.WF



Dimensional Drawing



Two-section Bobath table for rehabilitation and treatment of neuromotor disorders. Large surface with ample space for patient and therapist. Three face holes with plugs for maximum flexibility.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL02W.1.A66.WF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	190 × 120 cm
Height adjustment	49 – 93 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric height adjustment via perimeter foot switch
- Elegant design with rounded feet for added safety
- Backrest with 3 face holes and plugs
- Extra-wide surface: 120 cm
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 194 × 124 × 53 cm · Shipping weight approx. 101 kg · Shipping volume 1.28 m³