

EXAMINATION & TREATMENT COUCH · ELECTRIC

BROM PRO

Ref. EL03EPRO23SL.1.A66.GF



Three-section treatment table, ideal for medical centres, with face hole and plug. Electric height adjustment; back and leg sections adjustable via gas spring. Paper-roll holder as standard.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL03EPRO23SL.1.A66.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	3

DIMENSIONS

Lying surface (L × W)	195 × 70 cm
Height adjustment	45 – 95.5 cm
Backrest angle	+65° / +40°

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TIMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m³
Cover colours	blue or grey


CASTORS & SETUP

Castors	retractable, Ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Dimensional Drawing



Features & Equipment

- Electric height adjustment via foot switch
- Backrest with face hole and plug
- Elegant design with rounded feet for added safety
- Covered with very durable, fire-retardant vinyl

SHIPPING & PACKAGING Packaging size 186 × 70 × 54 cm · Shipping weight approx. 82 kg · Shipping volume 0.7 m³