

MOBILE STOOL WITH ARMREST

CORE

Ref. MA15AR.A66



Round stool with a 360° swivel armrest for easy movement and adaptation to various working positions. Movable swivel seat to relieve pressure on the legs – for comfort and ergonomics in continuous use.

PRODUCT

Manufacturer	COINFYCARE
Reference	MA15AR.A66
Category	Mobile stool with armrest
Version	Gas lift, swivel

DIMENSIONS

Seat diameter	40 cm
Height adjustment	48 – 64 cm
Base	Ø 58 cm

DRIVE & OPERATION

Height adjustment	gas lift
-------------------	----------


UPHOLSTERY & COVER

Padding thickness	7 cm
-------------------	------

PERFORMANCE & COMPLIANCE

Load capacity	160 kg
---------------	--------

Dimensional Drawing



Features & Equipment

- Round stool with 360° swivel armrest
- Movable swivel seat to relieve leg pressure
- For ergonomic continuous use

SHIPPING & PACKAGING Packaging size 57 × 25 × 55 cm · Shipping weight approx. 14 kg · Shipping volume 0.08 m³