

EXAMINATION & TREATMENT COUCH · ELECTRIC

ENID

Ref. EL02.1.A66.GF



Robust two-section table for intensive use in physiotherapy and medicine. Elegant design with a standard perimeter foot-switch system – ready to use immediately.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL02.1.A66.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	48 – 92 cm
Backrest angle	+70° / -20°

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	blue or grey


CASTORS & SETUP

Castors	retractable, Ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Dimensional Drawing



Features & Equipment

- Electric height adjustment via perimeter foot switch
- Elegant design with rounded feet for added safety
- Backrest with face hole and plug
- Very durable, fire-retardant vinyl upholstery

SHIPPING & PACKAGING Packaging size 199 × 74 × 51 cm · Shipping weight approx. 60 kg · Shipping volume 0.75 m³