

TREATMENT CHAIR

ENT

Ref. ELT35.4.A66.WF



Multifunctional electric treatment chair for diagnosis and therapy in ophthalmology, plastic surgery and ENT. Four motors control height, backrest and leg rest; incl. Trendelenburg position and presets for seated or flat position. Angle-adjustable headrest; high-resistance grey PVC upholstery.

PRODUCT

Manufacturer	COINFYCARE
Reference	ELT35.4.A66.WF
Category	Treatment chair
Version	Electric

DIMENSIONS

Overall length	198 cm
Width	90 cm
Height adjustment	59 – 89 cm

DRIVE & OPERATION

Height adjustment	electric
-------------------	----------

UPHOLSTERY & COVER

Padding thickness	75 mm
Foam density	36 kg/m ³



CASTORS & SETUP

Castors	ø 125 mm
---------	----------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year

Dimensional Drawing



Features & Equipment

- Electric Hi-Lo and angle operation via hand controller
- Luxurious, soft and comfortable head cushion
- Comfortable contoured table top with folding armrests
- Waste tray included
- Sleek design with rounded frame legs for added safety
- Padded, removable armrests
- Upholstery: PVC vinyl resistant to fire, stains, mould, oil, water and wear

SHIPPING & PACKAGING Packaging size 169 × 73 × 78 cm · Shipping weight approx. 146 kg · Shipping volume 0.96 m³