

EXAMINATION & TREATMENT COUCH · FIXED

FIX

Ref. FIXMT02.A66.WF



Fixed treatment couch – the workhorse of every practice. Flat-packed, assembled in two minutes; white frame with dark-grey upholstery. Hard-wearing, antibacterial upholstery, available in blue or grey.

PRODUCT

Manufacturer	COINFYCARE
Reference	FIXMT02.A66.WF
Category	Treatment / examination couch
Version	Fixed
Sections	2

DIMENSIONS


Lying surface (L × W)	195 × 70 cm
Height (fixed)	70 cm
Backrest angle	0° – 75°

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³

PERFORMANCE & COMPLIANCE

Load capacity	225 kg
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Dimensional Drawing

Features & Equipment

- Backrest angle: 0° to 75°

SHIPPING & PACKAGING Packaging size 197 × 73 × 14 cm · Shipping weight approx. 38 kg · Shipping volume 0.2 m³