

EXAMINATION & TREATMENT COUCH · ELECTRIC

GALLEY

Ref. EL05ARSL.2.A66.SH.GF



Seven-section physiotherapy table for versatile use. Adjustable armrests at the head end for the prone position. Two motors (height and electric kyphosis). Centre section with folding armrests (up to 90 cm wide) and paper-roll holder.

PRODUCT

Manufacturer	CONFYCARE
Reference	EL05ARSL.2.A66.SH.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	7

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	51 – 94 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TIMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	grey


CASTORS & SETUP

Castors	retractable, Ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Dimensional Drawing



Features & Equipment

- Electric height adjustment via perimeter control
- Elegant design with rounded feet for added safety
- Standard safety-lock system
- Armrests on the head section
- Folding armrests in the centre section (up to 90 cm wide)
- Covered with very durable, fire-retardant vinyl
- Upholstery colour: grey

SHIPPING & PACKAGING Packaging size 195 × 70 × 58 cm · Shipping weight approx. 76 kg · Shipping volume 0.79 m³