


EXAMINATION & TREATMENT COUCH · ELECTRIC

GUTH

Ref. EL022.1.A66.SH.GF



Dimensional Drawing



Two-section physiotherapy table with height-adjustable armrests at the head end for the ideal prone position. Headrest adjustable to positive and negative angles – for cervical and cranial techniques.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL022.1.A66.SH.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	48 – 92 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	blue or grey

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric height adjustment via perimeter foot switch
- Elegant design with rounded feet for added safety
- Head section with adjustable armrests
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 199 × 70 × 51 cm · Shipping weight approx. 64 kg · Shipping volume 0.71 m³