

EXAMINATION & TREATMENT COUCH · ELECTRIC

GUTH PRO

Ref. EL022PRO23SL.1.A66.SH.WF



Electric treatment table with 2 sections, armrests and a multi-position headrest for extended treatment and examination options. Ergonomic and safe for patient and practitioner.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL022PRO23SL.1.A66.SH.WF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	202 × 68 cm
Height adjustment	43 – 93 cm
Head section angle	–55° bis +45°

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER

Foam density	36 kg/m ³
--------------	----------------------


CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016

Dimensional Drawing



Features & Equipment

- Adjustable head section, angle –55° to +45°
- Electric Hi-Lo operation via perimeter foot bar
- Sleek design with rounded feet for added safety
- Head section with armrests for maximum comfort in the prone position
- Fully enclosed base covers all technical parts – easy to clean
- Upholstery: fire-retardant PVC vinyl per BS 5852, resistant to stains, mould and wear

SHIPPING & PACKAGING Packaging size 199 × 70 × 54 cm · Shipping weight approx. 82 kg · Shipping volume 0.75 m³