

EXAMINATION & TREATMENT COUCH · ELECTRIC

JULL PRO

Ref. EL03PRO23.1.A66.GF



Treatment table with a short head section as a solid, stable platform. Gas-assisted back and head sections and a quiet electric Hi-Lo lift system – a balanced ratio of performance, function and price.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL03PRO23.1.A66.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	3

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	45 – 95.5 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER

Foam density	36 kg/m ³
--------------	----------------------



CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
---------------	--------

Dimensional Drawing



Features & Equipment

- Adjustable angle: head section –55° to +45°, leg section 0° to +75°
- Electric Hi-Lo operation via perimeter foot bar
- Sleek design with rounded frame legs for added safety
- Head section with face hole and plug for maximum comfort in the prone position
- Fully enclosed base covers all technical components – easy to clean

SHIPPING & PACKAGING Packaging size 195 × 60 × 50 cm · Shipping weight approx. 94 kg · Shipping volume 0.58 m³