

EXAMINATION & TREATMENT COUCH · HYDRAULIC

KEND PRO

Ref. EL03PRO29.A66.WF



Table with a short head section as a solid, stable platform for professional treatments. Gas-assisted back and head sections and a quiet Hi-Lo lift system. Optimised access and improved ergonomics.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL03PRO29.A66.WF
Category	Treatment / examination couch
Version	Hydraulic (Hi-Lo)
Sections	3

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	45 – 95.5 cm

DRIVE & OPERATION

Height adjustment	hydraulic, side pedal
-------------------	-----------------------

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³



CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
---------------	--------

Dimensional Drawing



Features & Equipment

- Adjustable angle: head section –55° to +45°, leg section 0° to +75°
- Hydraulic Hi-Lo system via side foot pedal
- Sleek design with rounded frame legs for added safety
- Head section with face hole and plug for maximum comfort in the prone position
- Fully enclosed base covers all technical components – easy to clean

SHIPPING & PACKAGING Packaging size 195 × 66 × 48 cm · Shipping weight approx. 87 kg · Shipping volume 0.62 m³