


EXAMINATION & TREATMENT COUCH · HYDRAULIC

LING PRO

Ref. EL02PRO29.A66.GF



Dimensional Drawing



Robust two-section treatment table for demanding daily practice – whether in a busy physiotherapy clinic or a small general practice. Reliable, excellent performance.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL02PRO29.A66.GF
Category	Treatment / examination couch
Version	Hydraulic (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	45 – 95.5 cm
Backrest angle	+70° / –20°

DRIVE & OPERATION

Height adjustment	hydraulic, side pedal
Lifting time	12 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	blue or dark grey

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Hydraulic height adjustment of the table
- Sleek design with rounded feet for added safety
- Backrest with face hole and plug
- Very durable, fire-retardant vinyl upholstery

SHIPPING & PACKAGING Packaging size 199 × 70 × 52 cm · Shipping weight approx. 77 kg · Shipping volume 0.72 m³