

EXAMINATION & TREATMENT COUCH · HYDRAULIC

MAIT PRO

Ref. EL032PRO29.A66.GF



Very popular with physiotherapists: a quiet, durable hydraulic lift/lower system via side pedal and a wide range of motion of headrest and backrest. Optimal patient positioning with reduced strain on the practitioner.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL032PRO29.A66.GF
Category	Treatment / examination couch
Version	Hydraulic (Hi-Lo)
Sections	3

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	45 – 95.5 cm

DRIVE & OPERATION

Height adjustment	hydraulic, side pedal
-------------------	-----------------------

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³


CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
---------------	--------

Dimensional Drawing



Features & Equipment

- Adjustable angle: head section –55° to +45°, leg section 0° to +75°
- Hydraulic height adjustment (Hi-Lo) via side pedal
- Elegant design with rounded frame feet for added safety
- Head section with face hole and plug plus padded, height-adjustable armrests
- Fully enclosed base covers all technical parts – easy to clean
- Lift and lower time: 12 s

SHIPPING & PACKAGING Packaging size 199 × 70 × 50 cm · Shipping weight approx. 85 kg · Shipping volume 0.7 m³