

MASSAGE CHAIR

NOBU

Ref. CFMC03A.A66



Portable massage chair with ergonomic design for correct posture during treatment. Height and tilt of headrest, backrest, seat, armrests and knee supports adjustable in multiple ways. Foldable for transport and storage; easy to clean.


PRODUCT

Manufacturer	COINFYCARE
Reference	CFMC03A.A66
Category	Massage chair
Version	Foldable, manual

DIMENSIONS

Unfolded (W×D×H)	74 × 54 × 125 cm
Folded	26 × 135 cm

Dimensional Drawing



Features & Equipment

- Portable, foldable massage chair
- Height and tilt of headrest, backrest, seat, armrests and knee supports adjustable
- Ergonomic alignment of the user
- Easy to clean; simple transport and storage

SHIPPING & PACKAGING Packaging size 26 × 53 × 107 cm · Shipping weight approx. 13 kg · Shipping volume 0.15 m³