


EXAMINATION & TREATMENT COUCH · ELECTRIC

NOTT PRO

Ref. EL05PRO23SL.2.A66.SH.GF



Dimensional Drawing



Five-section physiotherapy table for versatile use. Height-adjustable armrests at the head end for the prone position. Two motors: height adjustment and electric kyphosis. Paper-roll holder as standard.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL05PRO23SL.2.A66.SH.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	5

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	45 – 95.5 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	blue or grey

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric height adjustment via perimeter foot switch
- Adjustable centre section (electric kyphosis)
- Elegant design with rounded feet for added safety
- Backrest with face hole and plug
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 193 × 70 × 51 cm · Shipping weight approx. 98 kg · Shipping volume 0.69 m³