

GYNAECOLOGICAL EXAMINATION CHAIR

OB

Ref. EL3603.2.A66.WF



Electric gynaecological examination chair for general examinations – indispensable in clinics, health centres and practices. Two motors for seat height and tilt ensure optimal handling; grey medical PVC upholstery, highly resistant and easy to clean.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL3603.2.A66.WF
Category	Gynaecological examination chair
Version	Electric

DIMENSIONS

Overall length	185 cm
Width	60 cm
Height adjustment	62 – 104 cm

DRIVE & OPERATION

Height adjustment	electric
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER


Padding thickness	5 cm
Foam density	36 kg/m ³

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year

Dimensional Drawing

Features & Equipment

- Electric seat height and tilt adjustment via pedal
- Elegant design with rounded legs for added safety
- Incl. foot rests and stirrups
- Incl. waste tray and side rails
- Filler section for general use
- Upholstery: high-resistance PVC vinyl, resistant to fire, stains, mould, oil, water and wear

SHIPPING & PACKAGING Packaging size 148 × 73 × 68 cm · Shipping weight approx. 89 kg · Shipping volume 0.73 m³