


EXAMINATION & TREATMENT COUCH · ELECTRIC

PEHR PRO

Ref. EL11PRO23SL.1.A66.WF



Dimensional Drawing



Two-section treatment table designed for physiotherapy. Headrest with face cushion and armrests beneath it, positive/negative adjustable. Contoured main section for easy patient access. Paper-roll holder as standard.

PRODUCT

Manufacturer	CONFY CARE
Reference	EL11PRO23SL.1.A66.WF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	197 × 70 cm
Height adjustment	48 – 92 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³
Cover colours	blue or grey

CASTORS & SETUP

Castors	retractable, Ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric height adjustment via perimeter foot switch
- Elegant design with rounded feet for added safety
- Tilttable headrest and armrests beneath it
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 193 × 70 × 50 cm · Shipping weight approx. 75 kg · Shipping volume 0.68 m³