

PODIATRY TREATMENT CHAIR

PODO

Ref. auf Anfrage



Electric podiatry treatment chair with three sections and padded armrests. Backrest and leg section widely adjustable for optimal treatment positions in podiatry and foot care. Electric height adjustment and comfortable, easy-care upholstery.

PRODUCT

Manufacturer	COINFYCARE
Reference	on request
Category	Podiatry treatment chair
Version	Electric


DIMENSIONS

Overall length	190 cm
Width (max.)	91 cm
Height adjustment	56 – 100 cm

DRIVE & OPERATION

Height adjustment	electric
-------------------	----------

Dimensional Drawing



Features & Equipment

- Three sections: back, seat and leg
- Adjustable backrest (up to 90°)
- Multi-position leg section (25° to -80°)
- Padded, folding armrests
- Head section with neck cushion
- Smooth-running castors
- Hard-wearing, easy-care PU upholstery
- Ideal for podiatry and foot care

SHIPPING & PACKAGING Packaging data on request