


EXAMINATION & TREATMENT COUCH · ELECTRIC

RILA

Ref. EL032.1.A66.SH.GF


Dimensional Drawing



Five-section couch for physiotherapy and osteopathy. Height-adjustable armrests at the head end, positive/negative headrest, fixed centre section. Paper-roll holder as standard.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL032.1.A66.SH.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	5

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	48 – 92 cm
Backrest angle	+80°

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³

CASTORS & SETUP

Castors	retractable, Ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Five-section couch for physiotherapy and osteopathy
- Head section with adjustable armrests
- Head section with face hole and plug; fixed centre section
- Fire-retardant, highly wear-resistant vinyl upholstery
- Paper-roll holder as standard

SHIPPING & PACKAGING Packaging size 194 × 70 × 53 cm · Shipping weight approx. 65 kg · Shipping volume 0.72 m³