

EXAMINATION & TREATMENT COUCH · ELECTRIC

ROTH

Ref. EL022AR.1.A66.GF



Two-section couch designed for osteopathy, with face hole and plug. Folding lateral armrests on the main section (+80 cm width). Headrest adjustable via gas spring, positive/negative. Paper-roll holder as standard.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL022AR.1.A66.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	2

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	48 – 92 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability
Lifting time	30 seconds

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	70 kg/m ³


CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year on frame, electrical components and upholstery (except in case of misuse)

Dimensional Drawing



Features & Equipment

- Electric height adjustment via perimeter control
- Elegant design with rounded feet for added safety
- Headrest with adjustable armrests
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 199 × 70 × 51 cm · Shipping weight approx. 64 kg · Shipping volume 0.71 m³