


CHIROPRACTIC COUCH · ELECTRIC

SCALL PRO

Ref. EL08BPRO23.1.A66.GF



Dimensional Drawing



Chiropractic table with drop system – versatile and economical. Cervical, thoracic, lumbar and pelvic drop sections, manual with tension control. Robust structure and electric lift system for safe, comfortable and effective treatments.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL08BPRO23.1.A66.GF
Category	Chiropractic manipulation table
Version	Electric (Hi-Lo)

DIMENSIONS

Overall length	185 – 200 cm
Width	65 / 55 cm
Height adjustment	56 – 90 cm

DRIVE & OPERATION

Height adjustment	electric (Hi-Lo), foot bar
Motor	TIMOTION, high reliability

UPHOLSTERY & COVER

Padding thickness	7.5 cm
Foam density	100 kg/m ³

CASTORS & SETUP

Levelling feet	yes, for uneven floors
----------------	------------------------

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	ISO 13485 / MDR-CE / FDA-registered
Warranty	2 years on frame, electrical components and upholstery (except in case of misuse)

Features & Equipment

- Electric Hi-Lo system via foot bar
- Head section with padded platform and adjustable width
- Head elevation via gas spring ($\pm 20^\circ$), 10 cm height adjustment
- Lateral head movement $\pm 15^\circ$
- Padded footrest, height-adjustable and extendable
- Four drop sections: cervical, thoracic, lumbar and pelvic
- 4 retractable castors with central locking for easy transport

SHIPPING & PACKAGING Packaging size 192 × 62 × 51 cm · Shipping weight approx. 92 kg · Shipping volume 0.61 m³