

CHIROPRACTIC COUCH · ELECTRIC

SIDO PRO

Ref. EL07PRO23.2.A66.WF



Versatile chiropractic manipulation table with maximum flexibility. For mobilisation, manipulation, massage and movement therapy; multi-angle adjustment for physiotherapy and rehabilitation. Exclusive functions such as tilt, inversion and lower lumbar rotation expand the therapeutic possibilities.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL07PRO23.2.A66.WF
Category	Chiropractic manipulation table
Version	Electric (Hi-Lo)


DIMENSIONS

Overall length	210 cm
Width	63 / 55 cm
Height adjustment	52 – 90 cm

DRIVE & OPERATION

Height adjustment	electric (Hi-Lo), foot bar
Motor	TiMOTION, high reliability

Dimensional Drawing



Features & Equipment

- Electric Hi-Lo system via foot bar; angle adjustment of upper body via button
- Armrests with dual-position adaptability
- Tilt function for assisted patient positioning
- Head section with width-adjustable padded platform
- Head elevation via gas spring (–15° to +35°), height adjustment 10 cm, lateral movement ±20°
- Two drop sections: cervical and thoracic; manual lumbar traction mechanism
- Vertical tension control and horizontal flexion of the lower body via T-bar
- Rotation of the lower body from 0° to 20°
- Dual-function extender as extension platform and step

SHIPPING & PACKAGING Packaging size 199 × 70 × 54 cm · Shipping volume 0.75 m³