


EXAMINATION & TREATMENT COUCH · ELECTRIC

STILL

Ref. EL052PRO23SL.2.A66.GF



Dimensional Drawing



Seven-section treatment table with powerful, quiet electric height adjustment (Hi-Lo) via perimeter foot bar and electric postural drainage via side button. Wide range of motion of the head and back sections.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL052PRO23SL.2.A66.GF
Category	Treatment / examination couch
Version	Electric (Hi-Lo)
Sections	7

DIMENSIONS

Lying surface (L × W)	195 × 66 cm
Height adjustment	51 – 94 cm

DRIVE & OPERATION

Height adjustment	electric, perimeter foot switch
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³

CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016

Features & Equipment

- Powerful, quiet electric height adjustment (Hi-Lo) via perimeter foot bar
- Electric postural drainage via side button
- Wide range of motion of head and back sections
- Seven sections for versatile use
- Fire-retardant, highly wear-resistant vinyl upholstery

SHIPPING & PACKAGING Packaging size 205 × 50 × 75 cm · Shipping weight approx. 101 kg · Shipping volume 0.77 m³