

BLOOD-DONATION CHAIR

XTRACT

Ref. EL35B.1.A66



Blood-donation couch for the donor's safety, comfort and calm. The human-centred design eases mental tension and makes the donation process simpler and more pleasant.

PRODUCT

Manufacturer	COINFYCARE
Reference	EL35B.1.A66
Category	Blood-donation chair
Version	Electric

DIMENSIONS

Overall length	190 cm
Width	60 cm
Height adjustment	56 – 100 cm

DRIVE & OPERATION

Height adjustment	electric
Motor	TiMOTION, high reliability

UPHOLSTERY & COVER

Padding thickness	5 cm
Foam density	36 kg/m ³


CASTORS & SETUP

Castors	retractable, ø 75 mm
Levelling feet	yes, for uneven floors

PERFORMANCE & COMPLIANCE

Load capacity	250 kg
Certification	CE, FDA, ISO 13485:2016
Warranty	1 year

Dimensional Drawing



Features & Equipment

- Electric Hi-Lo operation via perimeter bar
- Elegant design with rounded corners for added safety
- Comfortable flip-up armrests
- Removable padded headrest
- Upholstery: high-resistance PVC vinyl, resistant to fire, stains, mould, oil, water and wear

SHIPPING & PACKAGING Packaging size 174 × 73 × 63 cm · Shipping weight approx. 76 kg · Shipping volume 0.8 m³