




ROTHBAND



X-RAY GLASSES CATALOGUE


GLOBAL LEADERS IN RADIATION PROTECTION SOLUTIONS

LENS TYPE




PLAIN LENS

A non-prescription lens designed for eye protection and comfort. Ideal for users who do not require vision correction but still need reliable protection in clinical or radiation environments.




SINGLE VISION

A lens with a single, consistent prescription across the entire surface. Suitable for correcting either distance or near vision, depending on the user's requirements.




READER ONLY

A near-vision lens specifically designed for close-up tasks such as reading or detailed work. Ideal for users who require magnification for short working distances.



BIFOCAL

A dual-purpose lens featuring two distinct optical zones: one for distance vision and a lower segment for near tasks. Enables clear vision at two focal lengths within one lens.




VARIFOCAL

A multi-focal lens offering a smooth transition between distance, intermediate, and near vision. Provides seamless vision across all working distances without visible segment lines.

LENS COATINGS

Select only ONE, for Plain Lens, Single Vision, or Varifocal Prescriptions Only (Not available on Bifocal)



ANTI-REFLECTIVE

Recommended if working under fluorescent lights or with computer screens and devices




ANTI-FOG PERMANENT

Recommended if working with a mask or in a fogging environment.




LENS CLEANING / DISINFECTANT



UNCOTAED LENSES

- Clean with anti-bacterial detergent or hand soap with no microbeads or grit added.
- Gently rub lenses in straight lines under clear running water.
- Dry with soft, lint-free cloth.
- You can also use alcohol based lens wipes.



ANTI-REFLECTIVE LENSES

- **DO NOT USE ALCOHOL.**
It will degrade the anti-reflective coating.
- Clean with anti-bacterial detergent or hand soap with no microbeads or grit added.
- Gently rub lenses in straight lines under clear running water.
- Dry with soft, lint-free cloth.



ANTI-FOG LENSES

- **DO NOT USE ALCOHOL.**
It will degrade the anti-fog coating.
- Clean with anti-bacterial detergent or hand soap with no microbeads or grit added.
- You can also use non-alcohol based lens wipes.
- Use a light pressure in straight lines, thoroughly cleaning all surfaces of lenses and frame.
- Dry with microfiber cleaning cloth.



DON'TS

- **DO NOT AUTO-CLAVE.**
- **NEVER USE AMMONIA OR CHLORINE BASED CLEANERS.**
- **DO NOT USE CIRCULAR MOTIONS WHILE CLEANING, AS IT CAN DAMAGE APPLIED COATINGS.**




NIKE BRAZEN

SKU: GLNIKESS



FRAME INFORMATION



The Nike Brazen features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. The wraparound-style frame is available with prescription and includes rubberised nose pads and temple bars, plus side lead protection (also available without side lead protection).

Weight: 80g

Available in: Black or matte black with military blue

SPECIFICATION

Size: 60-18-134

Weight: 80g

Protection level: 0.75mm Pb LE
0.50mm Pb side

Colour: Black or matte black with military blue

Material: Plastic

Shape: Wrap Around

Prescription Available: Yes



WSR Medical Solutions Ltd T/A Rothband

Unit 2 Hargreaves Mill, Hargreaves Street, Haslingden, Lancashire, BB4 5RQ

T

+44 (0) 1706 830086

W


www.rothband.com

OAKLEY GASCAN

SKU: GLOGC



FRAME INFORMATION



The Oakley Gascan features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This wraparound-style frame is available with prescription and features a single bridge. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 80g

Available in: Polished black or matte black

SPECIFICATION

Size:	60-15-127
Weight:	80g
Protection level:	0.75mm Pb LE
Colour:	Polished black or matte black
Material:	Plastic
Shape:	Wrap Around
Prescription Available:	Yes




MODEL 33

FIT-OVER FRAME

SKU: GLFT-M



FRAME INFORMATION



The Medium Fit-Over features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This non-prescription frame features adjustable temples and integrated side shields.

Weight: 90g

Available in: Black or tortoise

SPECIFICATION

Size:	60-16-127
Weight:	90g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black or tortoise
Material:	Plastic
Shape:	Rectangular / Fit Over
Prescription Available:	No




MODEL 38

FIT-OVER FRAME

SKU: GLFT-L



FRAME INFORMATION



The Large Fit-Over features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This non-prescription frame features adjustable temples and integrated side shields.

Weight: 90g

Available in: Black or tortoise

SPECIFICATION

Size:	65-20-125
Weight:	90g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black or tortoise
Material:	Plastic
Shape:	Rectangular / Fit Over
Prescription Available:	No




MODEL 52

SKU: GLWR



FRAME INFORMATION



The GLWR features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features adjustable nose pads and integrated side shields.

Weight: 90g

Available in: Black, tortoise, or blue

SPECIFICATION

Size:	53-15-120
Weight:	90g
Protection level:	0.75mm Pb front 0.75mm Pb side
Colour:	Black, tortoise or blue
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL 100

SKU: GL100



FRAME INFORMATION



The GL100 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features adjustable nose pads, leaded side shields, and spring hinges.

Weight: 89g

Available in: Black or brown

SPECIFICATION

Size: 56-21-145 | 52-19-140

Weight: 83g | 89g

Protection level: 0.75mm Pb front
0.50mm Pb side

Colour: Black or brown

Material: Metal

Shape: Aviator

Prescription Available: Yes




MODEL 206

SKU: GL206



FRAME INFORMATION



The GL206 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads, rubberised temples, integrated side shields, and a removable dust dam.

Weight: 75g

Available in: Orange with brown with clear side shields, yellow with black with clear side shields, or yellow with black with smoke grey side shields

SPECIFICATION

Size:	50-19-117
Weight:	75g
Protection level:	0.75mm Pb front 0.5mm Pb side
Colour:	Orange with brown with clear side shields, yellow with black with clear side shields, or yellow with black with smoke grey side shields
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL 208

SKU: GL208



FRAME INFORMATION



The GL208 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses have clear CE-certified lenses. This frame is available with prescription and features rubberised temple bars. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 72g

Available in: Black, blue with yellow, pink or silver

SPECIFICATION

Size: 59-20-130

Weight: 72g

Protection level: 0.75mm Pb

Colour: Black, blue with yellow, pink or silver

Material: Plastic

Shape: Wrap Around

Prescription Available: Yes




MODEL 400

SKU: GL400



FRAME INFORMATION



The GL400 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features silicone nose pads, spring hinges, and removable side shields.

Weight: 64g

Available in: Gunmetal or black

SPECIFICATION

Size:	54-19-145
Weight:	64g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Gunmetal or black
Material:	Metal
Shape:	Rectangular
Prescription Available:	Yes




MODEL 1171

SKU: RG1171



FRAME INFORMATION



The RG1171 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This non-prescription frame features a foam eyecup, rubberised nose pads, and rubberised temple tips. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 69g

Available in: Red or black

SPECIFICATION

Size:	58-15-125
Weight:	69g
Protection level:	0.75mm Pb
Colour:	Red or black
Material:	Plastic
Shape:	Wrap Around
Prescription Available:	No




MODEL 1205

SKU: GL1205



FRAME INFORMATION



The GL1205 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This non-prescription frame features rubberised nose pads and temple bars. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 71g

Available in: Blue or black

SPECIFICATION

Size:	64-16-140
Weight:	71g
Protection level:	0.75mm Pb
Colour:	Blue or black
Material:	Plastic
Shape:	Wrap Around
Prescription Available:	No




MODEL 9388

SKU: GL9388



FRAME INFORMATION



The GL9388 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features adjustable temples and permanent side shields.

Weight: 73g

Available in: Black

SPECIFICATION

Size:	63-14-120
Weight:	73g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black
Material:	Plastic
Shape:	Oversized
Prescription Available:	Yes




MODEL 15011

SKU: GL15011



FRAME INFORMATION



The GL15011 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads, rubberised temple bars, adjustable nose pads, and permanent side shields.

Weight: 88g

Available in: Black with clear, black with green, crystal red, crystal clear, navy with clear, crystal purple, cyan blue, or crystal faded grey

SPECIFICATION

Size:	57-15-125
Weight:	88g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black with clear, black with green, crystal red, crystal clear, navy with clear, crystal purple, cyan blue, or crystal faded grey
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL 17012

SKU: GL17012



FRAME INFORMATION



The GL17012 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubber temple bars and side shields.

Weight: 72g

Available in: Black

SPECIFICATION

Size:	50-20-135
Weight:	72g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL 20020

SKU: GL20020



FRAME INFORMATION



The GL20020 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised temple bars and integrated side shields.

Weight: 81g

Available in: Clear blue

SPECIFICATION

Size:	52-18-128
Weight:	81g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Clear blue
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL EGM

SKU: GLEGM



FRAME INFORMATION



The GLEGM features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads and rubberised temples. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 88g

Available in: Black, grey, or red

SPECIFICATION

Size:	53-16-122
Weight:	88g
Protection level:	0.75mm Pb
Colour:	Black, grey or red
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL EX601

SKU: EX601



FRAME INFORMATION



The EX601 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features a foam eyecup, rubberised nose pads, rubberised temples, and leaded side shields.

Weight: 77g (89g with foam gasket)

Available in: Black, black with foam gasket, or matte apple green with foam gasket

SPECIFICATION

Size:	54-17-130
Weight:	77g (89g with foam gasket)
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black, black with foam gasket, or matte apple green with foam gasket
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL F10

SKU: F10-EG



FRAME INFORMATION



The F10 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads, rubberised temple bars, and integrated side shields.

Weight: 93g

Available in: Clear, black with blue, or black with yellow

SPECIFICATION

Size:	55-19-125
Weight:	93g
Protection level:	0.75mm Pb front 0.50mm Pb side

Colour:	Clear, black with blue, or black with yellow
---------	--

Material:	Plastic
-----------	---------

Shape:	Rectangular
--------	-------------

Prescription Available:	Yes
-------------------------	-----




MODEL P820

SKU: GLP820



FRAME INFORMATION



The GLP820 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This non-prescription frame features rubberised nose pads. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 78g

Available in: Black or silver

SPECIFICATION

Size:	63-20-127
Weight:	78g
Protection level:	0.75mm Pb
Colour:	Black or silver
Material:	Plastic
Shape:	Wrap Around
Prescription Available:	No




MODEL RG12-PTG-038

FIT-OVER FRAME

SKU: RG12-PTG-038



FRAME INFORMATION



The RG12 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This non-prescription fit-over frame features a flip lens made of lead acrylic, which can be easily flipped up when radiation protection is not needed. These glasses fit comfortably over most prescription eyewear and offer 0.12mm Pb protection in the flip lens.

Weight: 70g

Available in: Black

SPECIFICATION

Size:	150-18-125
Weight:	70g
Protection level:	0.12mm Pb in the flip
Colour:	Black
Material:	Plastic
Shape:	Fit Over
Prescription Available:	No




MODEL Q100

SKU: Q100



FRAME INFORMATION



The Q100 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads and temple bars.

Weight: 69g

Available in: Black

SPECIFICATION

Size:	50-20-119
Weight:	69g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black
Material:	Plastic
Shape:	Wrap Around
Prescription Available:	Yes




MODEL Q200

SKU: Q200



FRAME INFORMATION



The Q200 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads and temple bars.

Weight: 74g

Available in: Black, purple, or matte blue with pink

SPECIFICATION

Size:	54-18-127
Weight:	74g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black, purple, or matte blue with pink
Material:	Plastic
Shape:	Rectangular, Wrap Around
Prescription Available:	Yes




MODEL Q300

SKU: Q300



FRAME INFORMATION



The Q300 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features rubberised nose pads, rubberised temple bars, and side shields.

Weight: 69g

Available in: Black

SPECIFICATION

Size:	55-17-120
Weight:	69g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black
Material:	Plastic
Shape:	Rectangular, Wrap Around
Prescription Available:	Yes




MODEL RK2

SKU: GLRK2



FRAME INFORMATION



The GLRK2 goggles feature high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This frame is available with prescription and features a foam eyecup and strap.

Weight: 85g

Available in: Black

SPECIFICATION

Size: 60-16-STRAP

Weight: 85g

Protection level: 0.75mm Pb

Colour: Black

Material: Plastic

Shape: Goggles

Prescription Available: Yes



WSR Medical Solutions Ltd T/A Rothband

Unit 2 Hargreaves Mill, Hargreaves Street, Haslingden, Lancashire, BB4 5RQ

T

+44 (0) 1706 830086

W


www.rothband.com

MODEL T9559

SKU: GLT9559



FRAME INFORMATION



The GLT9559 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This frame is available with prescription and features silicone nose pads, rubberised temples, removable foam support, and integrated side shields.

Weight: 76g

Available in: Clear neon green, clear teal, clear pink, or clear black

SPECIFICATION

Size:	58-17-128
Weight:	76g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Clear neon green, clear teal, clear pink, or clear black
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes




MODEL TP280

SKU: TP280



FRAME INFORMATION



The TP280 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection and clear CE-certified lenses. This non-prescription frame features a rubberised nose pad and temple bars. Side lead is not required, as the lens curvature provides 0.75mm Pb lateral protection.

Weight: 87g

Available in: Black with camo

SPECIFICATION

Size: 65-17-127

Weight: 87g

Protection level: 0.75mm Pb

Colour: Black with camo

Material: Plastic

Shape: Wrap Around

Prescription Available: No




MODEL X25

SKU: GLX25



FRAME INFORMATION



The GLX25 features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. These glasses also provide 0.50mm Pb lateral protection and have clear CE-certified lenses. This non-prescription frame features leaded side shields.

Weight: 94g

Available in: Black or tortoise

SPECIFICATION

Size:	56-15-135
Weight:	94g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black or tortoise
Material:	Plastic
Shape:	Square
Prescription Available:	No




MODEL T9538S

SKU: GLT9538S



FRAME INFORMATION



The GLT9538S features high-quality, distortion-free SF-6 Schott glass radiation-reducing lenses with 0.75mm Pb lead equivalency protection. This frame is available with prescription and features removable 0.50mm Pb side shields and clear CE-certified lenses.

Weight: 76g

Available in: Black, clear, or brown

SPECIFICATION

Size:	55-19-134
Weight:	76g
Protection level:	0.75mm Pb front 0.50mm Pb side
Colour:	Black, clear, or brown
Material:	Plastic
Shape:	Rectangular
Prescription Available:	Yes



LEAD SIDE SHIELD

UNIVERSAL RADIATION PROTECTIVE SIDE SHIELDS

SKU: GL-SS



The GL-SS side shields provide innovative lateral eye radiation protection. They are made from nylon-covered leaded vinyl and offer excellent hygiene properties. These side shields can be easily added to any eyewear and replaced when necessary. They provide 0.25mm Pb lead equivalency.

Available in: Black, purple, or blue

SPECIFICATION

Size:	2" H x 3.75" W
Weight:	50g - pair
Protection level:	0.25 mm Pb LE
Colour:	Black, purple, or blue
Material:	Leaded Vinyl
Features:	Easily added to any eyewear, optimum hygiene



TERMS and Conditions

1. Parties

- a) The "Company" means WSR Medical Solutions Ltd.
- b) The "Customer" means the person, firm or company buying goods from the Company to whom the goods are to be invoiced (unless the context otherwise requires) and any person to whom or for the benefit of whom the goods may be delivered.
- c) The "Goods" means the goods described on our invoice.

2. Application Of Conditions

- a) No contract shall exist between the Company and the Customer except upon the conditions set out herein and shall not be varied or excluded by the terms and conditions of any Order or any other terms and conditions sought to be imposed by the Customer including but without limitation to the Customer's Conditions of Purchase. These Terms and Conditions embody the entire understanding of the parties and supersede all and any prior promises, representations, undertakings or implications. No variation of these Terms and Conditions is permitted unless expressly accepted by a director of the Company in writing.
- b) If the Customer shall make default in or commit a breach of any of its obligations under this contract or any other Order or Contract with the Company, the Company shall at any time thereafter be entitled to determine the contract and without prejudice to its other rights hereunder recover from the Customer any loss arising.
- c) No statement, illustration or drawing in any circular, catalogue, advertisement or other communication in whatsoever form shall be deemed to imply any warranty or condition unless the same is specifically incorporated in this contract.

3. Orders

- a) Orders with a value of less than £20 (excluding VAT, net discounts) will attract a £2.50 low order surcharge to cover administration costs.
- b) Orders requesting the Company to provide the goods may be submitted or varied either orally or in writing, to either the Sales Representative, Customer Care department or other authorised Sales/Marketing personnel. Any orders submitted or variations made orally must be confirmed in writing within seven days.
- c) All orders to be forwarded to the address below.

4. Quotations

Any estimate or quotation shall comprise an invitation to treat and will remain open for a period of thirty days only from the date thereof, unless varied by the Company in writing. Any order issued by the Customer is subject to acceptance by the Company and a contract will only be formed when the Company has accepted the Customers offer to buy either in writing or by the processing of the order.

5. Specification

The Goods are supplied in accordance with specifications (if any) submitted to the Customer and any additions and alterations thereto shall be the subject of an extra charge. The Company reserves the right to change the construction or design of the Goods or the specification thereof and the Customer shall accept the Goods so changed in fulfillment of the Order.

6. Prices

- a) All prices quoted are exclusive of Value Added Tax which shall be payable by the Customer at the rate prevailing at the date of the applicable tax point.
- b) The Company reserves the right, by written notice at any time before despatch, to vary the price of the Goods to take account of any increases in the cost of raw materials, manufacture, packing, transport or wages or otherwise.
- c) Prices charged will be those effective at the date of Order as per the current price list.

7. Payment

- a) Unless the Company otherwise agrees, payment is due in full within 30 days net of the invoice date to WSR Medical Solutions Ltd., at the address below.
- b) If payment is not made in accordance with paragraph (a), or if at any time the credit standing of the Customer has (in the opinion of the Company) been impaired, the Company may refuse delivery of any Goods required under the Contract until arrangements as to payment or credit have been established which are satisfactory to the Company.
- c) Time of payment shall be of the essence and if payment is not made on the due date, the Company reserves the right to charge interest on all overdue payments and outstanding balances at the rate (both before and after judgement) of 5% above the Bank of England base rate in force on the due date for payment (under clause 7 a)) such interest being deemed to accrue on a daily basis from the due date for payment until receipt of payment by the Company.
- d) The Customer shall not be entitled to withhold payment either wholly or in part once the same has become due except where expressly agreed with the Company.

8. Delivery / Carriage Charges

- a) Unless the Company otherwise agrees, delivery shall be made to the delivery point specified in the Order.
- b) Time of delivery is not of the essence. Stated delivery times are no more than an estimate on the part of the Company and shall not be binding upon the Company which shall not be liable for any costs or damage caused by reason of any delay in delivery.
- c) The Company reserves the right to make delivery by instalments and to tender a separate invoice in respect of each instalment.
- d) The Company will usually make a charge for Carriage & packing and unless otherwise instructed by the Customer will choose the method of shipment.

9. Cancellations

- a) No notice of cancellation of an Order shall be valid unless given in writing and agreed to by an authorised representative of the Company.
- b) If written notice of cancellation is received by the Company, the Company reserves the right to levy a cancellation charge up to 25% of the amount of the Order which shall be paid by the Customer to the Company forthwith.
- c) If the Contract provides for delivery by instalments, each instalment shall be deemed to be the subject of a separate contract and non delivery or delay in the delivery of any instalment shall not entitle the Customer to repudiate the contract as regards the balance thereof.

10. Passing Of Title

- a) The legal title to the Goods shall not pass to the Customer until all sums due or payable by the Customer to the Company whether in respect of the Goods or otherwise howsoever shall have been received by the Company. Until such time the Customer shall hold the Goods as bailee of and in a fiduciary capacity for the Company and shall further:
 - i) store the Goods separately from any other goods so as to be identifiable as the property of the Company;
 - ii) be responsible for any loss or damage to the Goods howsoever caused; iii) except as provided in sub-clause (b) hereof, not sell, pledge, assign, charge or otherwise dispose of the Goods or any interest therein;
 - iv) deliver up the Goods to the Company at any time on demand and/or permit the Company, its servants or agents, at any time without notice to enter upon the Customer's premises and retake possession of the Goods.
- b) The Customer shall be entitled to sell in the ordinary course of its business any of the Goods, in which title shall not have passed to it, on condition that the Customer shall hold on trust for and (on demand made by the Company) pay or transfer to the Company the proceeds of such sale and all claims that the Customer may have against its purchaser as the result of such sale.
- c) The Customer (if so requested by the Company) and any liquidator or receiver of the assets of the Customer shall pay into a separate bank account for the sole benefit of the Company all sums received in respect of any sale under the powers out in sub-clause (b) above, by the Customer or the liquidator or receiver of any of the Goods in which title has not passed to the Customer.

11. Passing Of Risk

The risk in the Goods shall pass to the Customer when the Goods are tendered for delivery at the delivery point specified in the Order. In the event that the Goods are collected by the Customer from the Customers works the risk of the Goods shall pass to the Customer when they are loaded onto the vehicle of the Customer or of his carrier or other agent.

12. Warranty

- a) The Company warrants that the Goods shall at the time of delivery correspond to the specifications published by the Company when used for purposes for which goods of that type are ordinarily used. In the absence of any such published specifications applicable to the Goods, the Company warrants that the Goods shall at the time of delivery be of merchantable quality. If any Goods do not conform to that warranty the Company will at its option:
 - i) replace the Goods found not to conform to the warranty;
 - ii) bring the Goods into conformity with the published specifications of the Company or (as the case may be) merchantable quality; or
 - iii) take back the Goods found not to conform to the warranty and refund the appropriate part of the purchase price.

Provided that:

- aa) the liability of the Company shall in no event exceed the purchase price of the Goods; and
- bb) performance of any one of the above options (as limited by (iii) (aa) above) shall constitute an entire discharge of the Company's liability under this warranty.

b) The foregoing warranty is conditional upon:

- i) the Customer giving written notice to the Company of the alleged defect in the Goods such notice to be received by the Company within seven days or the time when the Customer discovers or ought to have discovered the defect;
- ii) the Customer affording the Company a reasonable opportunity to inspect the Goods;
- iii) the Customer making no further use of the Goods that are alleged to be defective, or any plant machinery or other equipment into which the same shall have been incorporated after the time at which the Customer discovers or ought to have discovered that they are defective.

- c) No representation or warranty is given as to the suitability or fitness of the Goods for any particular purpose and the Customer shall satisfy himself in this respect and shall be totally responsible therefore.
- d) Save as provided in paragraph (a) hereof and in Section 12 of the Sales of Goods Act, 1979:
 - i) all conditions and warranties express or implied are to the extent permitted by law hereby expressly excluded;
 - ii) the Company shall be under no liability for any loss or damage (other than in respect of personal injury or death) whether direct, indirect or consequential howsoever arising which may be suffered by the Customer.
 - e) Without prejudice to the generality of paragraph (c) all recommendations and advice given by or on behalf of the Company to the Customer as to the methods of storing, applying or using the Goods, the purposes to which may be applied and the suitability of using the Goods in any manufacturing process or in conjunction with any other materials are given without liability on the part of the Company its servants or agents and where the Goods are intended for use in connection with any foods, drugs, or other goods or substances whatever, or for the packaging or securing of goods or articles of any kind the Customer shall be responsible for ensuring that such articles will not be, or not liable to be adversely affected by contact with or proximity to the Goods supplied and all liability to the Customer or any third party in respect of any damage so caused is hereby expressly excluded.
 - f) In the event that, notwithstanding the provisions of paragraphs (c) and (d) above, the Company is found liable for any loss or damage suffered by the Customer, that liability shall in no event exceed the purchase price of the Goods.
 - g) To the extent that any provision herein contained or becomes unlawful by virtue of any statutory provisions or any other rule of law the words or phrases which so offend shall be deemed to have been deleted herefrom to the extent necessary to remove the invalidity, but no further or otherwise, and the remaining words, phrases and clauses shall remain operative, notwithstanding the deletion, and shall be construed, where possible, to give effect to the intentions of the parties as evidenced by the terms hereof.
 - h) No representation shall be taken to have been given or implied from anything said or written in the negotiations between the parties hereto or their representatives prior to any agreement entered into between the Company and the Customer but nevertheless if any representation shall subsequently be found or deemed to have been given it is hereby expressly excluded or restricted to the extent permitted by Section 3 of the Misrepresentation Act 1967 as amended by the Unfair Contract Terms Act 1977.

13. Force Majeure

The Company shall not be liable for any loss or damage caused by the delay in the performance or non-performance of any of its obligations hereunder where the same is occasioned by any cause whatsoever that is beyond the Company's control including but not limited to an act of god, war, civil disturbance, requisitioning, governmental or parliamentary restrictions, prohibitions or enactment's of any kind; import or export regulations, strike, lock-out or trade dispute (whether involving its own employees or those of any other person), difficulties in obtaining workmen or materials; breakdown of machinery, fire or accident. Should any such event occur the Company may cancel or suspend this contract without incurring any liability for any loss or damage thereby occasioned.

14. Notice Of Damage, Non-Delivery Or Wrong Delivery

- a) The Company shall not be liable to the Customer in the event that the Goods, at the time of delivery, do not comply with the warranty in Clause 12(a) and the non-compliance is one which should have been apparent on a careful visual inspection of the Goods unless the Customer shall have notified the Company and the Carrier within three days of delivery that such inspection was impossible and shall have confirmed such notification in writing received by the Company within seven days of delivery.
- b) The Company gives no guarantee of the exact quantity of goods delivered and in the event that quantity of Goods delivered does not correspond with the quantity stated on the Order the Customer shall only be liable to pay for the quantity delivered in the case of short-delivery and for the contract price in the event of over-delivery (subject in the latter case to permitting the Company to collect the surplus goods) provided that in the case of short delivery the Company is given an opportunity to inspect Goods before any use is made thereof or any alteration or modification is made thereto and that in no event shall such short or over-delivery entitle the Customer to damages or give the Customer the right to rescind.
- c) Claims for non delivery must be notified to the Company and the Carrier within twenty one days from date of invoice, failing which the Company shall be under no liability for non delivery.
- d) Any surplus goods delivered shall remain the property of the Company and the Customer shall take all reasonable precautions for the safe custody and protection of such surplus goods until the time of their removal by the Company.
- e) In the event of goods being wrongly delivered in no circumstances shall any of the Goods be returned to the Company without its prior written consent, otherwise liability is not accepted. Where the Company has made the error, full credit will be given on uplifted goods in good condition.
- f) The Company is under no obligation to accept goods for return, e.g. in the case of Customer ordering error, surplus stock, etc., but where it agrees to do so in writing a 15% handling charge will apply. Goods held by the Customer for more than 3 months from delivery will not normally be accepted for return. Only complete units of sale may be returned and all goods are to be in good condition. Any goods returned in unsaleable condition will be paid for in full by the Customer. The Company will, in all cases, organise the collection of goods.
- g) Refused Deliveries. The Company reserves the right to make a carriage charge plus 15% handling charge on all consignments tendered for delivery, which are refused without justification. All goods to be in good condition. Goods returned in an unsaleable condition will be paid for in full by the Customer. Goods should not be returned without the prior agreement of the Company.

15. Disputes And Set-Off

Any liability of the Company under the contract shall be subject to and conditional upon the due performance and observance by the Customer of all its obligations under these conditions and subject to these conditions the Customer shall not be entitled to withhold or delay payment or exercise any right or set-off whatsoever and howsoever arising or arisen which might otherwise be available to it.

16. Copyright

All copyrights shall remain the property of the Company alone unless otherwise expressly agreed in writing with the Customer.

17. Patent Rights

The Customer warrants that any instructions furnished or given by the Customer shall not be such as will cause the Company to infringe any letters patent or trade mark or other intellectual property right in execution of the Customer's order and shall indemnify the Company against all actions, processings, costs, claims and such demands arising therefrom to the extent that such instructions lead to such infringements.

18. Government Charges

The Customer will reimburse the Company for any increase in tax or Governmental charge, or for any new tax or Governmental charge hereafter becoming effective which the Company may be required to pay to any Government upon the sale, production, or transportation of the Goods and which has the effect of increasing the cost of the Goods sold hereunder.

19. Assignment

This contract is personal to the Customer and may only be assigned by it with the prior written consent of the Company.

20. Bankruptcy / Liquidation

If the Customer being an individual (or, when the customer is a firm, any partner in that firm) shall at any time become bankrupt, or having a receiving order made against them or shall make any composition or arrangement with or for the benefit of the creditors or, if the Customer being a company, an order for winding up is made or a resolution for winding up is passed or if a receiver or manager is appointed on behalf of a creditor or otherwise, the Company shall have the right to stop Goods in Transit, to suspend further deliveries and to determine any subsisting contract with the Customer and (in any case where title to the Goods has passed to the Customer but payment is due and has not been made) the Company shall have a lien over the Goods and right of resale provided that any such action shall be without prejudice to any right or claim the Company may have against the Customer.

21. Waiver

No neglect, delay or indulgence on the part of the Company in enforcing the terms and conditions of the contract shall prejudice the rights of the Company hereunder or be construed as a waiver thereof.

22. Terms Of Contract And Controlling Law

These conditions shall constitute the whole terms of any contract entered into between the Company and the Customer and shall be read and construed in accordance with the laws of England and all disputes arising in connection with the contract shall be submitted to the jurisdiction of the English courts.

23. Company Details

WSR Medical Solutions Ltd., Unit 2 Hargreaves mill, Hargreaves Street, Haslingden, Lancashire BB4 5RQ.
Tel: +44 (0) 1706 830086 - Email: sales@rothband.co.uk



WSR Medical Solutions Limited

Unit 2 Hargreaves Mill, Hargreaves Street,
Haslingden, Lancashire BB4 5RQ

T. +44 (0) 1706 830086 www.rothband.com
E. sales@rothband.com www.xray4vets.com

FOR TERMS AND CONDITIONS OF SALE, PLEASE VISIT: <https://www.rothband.com/terms-conditions-of-sale>

All of the product photos, colours and designs contained within this catalogue serve only as a representation. Colour and design deviations may occur as a result of the printing process and product development activities.

To the best of our knowledge, this catalogue has been put together accurately. Errors that occur as a result of product development, cannot be excluded.
All trademarks, logos and product names belong to their respective owners.

©2024 WSR Medical Solutions Limited

