

Table-Mounted Radiation Protection.

Specifications & Ordering Information




Table-Mounted Radiation Protection · 312/E-021



Contents

About This Product Range	3
Design & Key Features	4
Shield Configurations · 78–162 cm	5
Collision-Protected Shields · 92–162 cm	11
Top-Up Shields · Options, Heights & Compatibility	15
Head-End Shields	19
Foot-End Shield	23
Storage · Wall Brackets & Mobile Stands	24
Contact	27



About This Product Range


KENEX Lower-Body Shields considerably reduce the scattered radiation emitted by the table or the patient's body. They consist of a lightweight, durable structure with a reinforced double friction joint or single friction joint and a swivel arm carrying a flexible protective curtain with 0.5 mm Pb lead equivalent. Optional, semi-flexible top-up shields provide additional protection for the upper body.

- All shields can be mounted on either side of the table, so that the swivel arm always points towards the head end.
- The mounting bracket is pushed onto the table rail and quickly adapts to 25 mm or 28.6 mm high rails (max. 11 mm thick).
- Automatic wear compensation ensures constant joint friction – even with the table tilted in any direction.
- The hanging curtain forms a continuous protective zone for radiologists and cardiologists.
- The curtain wraps around the under-table X-ray tube without impeding movement.
- In combination with ceiling-mounted KENEX systems, upper- and lower-body protection complement each other to form an effective scattered-radiation barrier.




312/DS-039/6/3 · Lower-Body Shield with top-up shields

Design & Key Features



Integrated Accessory Rail

All shields from 100 cm width are supplied with two integrated accessory rails (10×25 mm mounted, 9.5×28.6 mm supplied). Control devices can be attached close to the user.




Stable Joint Pins

Internal stainless steel pins (10 mm Ø) hold the components securely and provide stable, permanently load-bearing pivot points.



Both-Side Mounting

All models can be quickly attached to the left or right of the table and moved to the other side if required.



Self-Adjusting Joint

A maintenance-free mechanism maintains the set friction and prevents unwanted curtain movement – for a long service life.



Rail-Width Adjustment

Secure attachment to 25 mm or 28.6 mm high table rails, configurable on site (max. 11 mm thick).



Curtain Construction

Flexible inner core with a robust outer shell (0.5 mm Pb) – hard-wearing, easy to clean, does not delaminate. Optionally as three overlapping panels (add-on /P).



Collision Protection

Patented mechanism: the swivel arm deflects upwards on collision with the C-arm – protecting shield and equipment.



Handle with Protective Collar

Handle made of impact-resistant, glass-fibre-reinforced plastic, protected by a solid collar around the handle shaft.

> Table-Mounted Shields Shield Configurations

The KENEX range for lower-body protection spans from 78 cm to 162 cm in width. Each shield is supplied with a flexible protective curtain with 0.5 mm Pb lead equivalent. Top-Up Shields extend the protection and, depending on the product, are already included or available separately (compatibility see page 16).

78 cm wide protective shields

Consisting of:

- Curtain Panel 78 cm wide × 70 cm high (optionally as three overlapping panels).
- Single Friction Joint for precise positioning.
- Top-Up Shields included with some products – compatibility see page 16.



78 cm wide Product Code	Curtain Panel		Top-Up Shields · 57 cm	
	single-piece	three overl.	Lift-off 312/DS/3,32	Fold-down 312/DS/3,33
312/E-005	●			
312/E-005/P		●		
312/E-006	●		●	
312/E-006/P		●	●	
312/E-007	●			●
312/E-007/P		●		●

● included · Further top-up shields available – compatibility see page 16

92 cm wide protective shields

Consisting of:


- Curtain Panel 92 cm wide × 70 cm high (optionally as three overlapping panels).
- Double Friction Joint for precise positioning.
- Top-Up Shields included with some products – compatibility see page 16.



312/E-011



312/E-012/P



312/E-010

92 cm wide Product Code	Curtain Panel		Top-Up Shields · 57 cm	
	single-piece	three overl.	Lift-off 312/DS/3.32	Fold-down 312/DS/3.33
312/E-010	●			
312/E-010/P		●		
312/E-011	●		●	
312/E-011/P		●	●	
312/E-012	●			●
312/E-012/P		●		●

● included · Further top-up shields available – compatibility see page 16

100 cm wide protective shields

Consisting of:


- Curtain Panel 100 cm wide × 70 cm high (optionally as three overlapping panels).
- Single Friction Joint for precise positioning.
- Two integrated Accessory Rails 35 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied); optionally 7.7×25 mm.
- Top-Up Shields included with some products – compatibility see page 16.



312/E-016



312/E-016/P



312/E-015

100 cm wide Product Code	Curtain Panel		Accessory Rail · 35 cm		Top-Up Shields · 57 cm	
	single-piece	three overl.	10×25 mm 9,5×28,6 mm	7.7×25 mm	Lift-off 312/DS/3,32	Fold-down 312/DS/3,33
312/E-015	●		●			
312/E-015/P		●	●			
312/E-015/8	●			●		
312/E-015/8P		●		●		
312/E-016	●		●		●	
312/E-016/P		●	●		●	
312/E-016/8	●			●	●	
312/E-016/8P		●		●	●	
312/E-017	●		●			●
312/E-017/P		●	●			●
312/E-017/8	●			●		●
312/E-017/8P		●		●		●

● included · Further top-up shields available – compatibility see page 16

114 cm wide protective shields

Consisting of:


- Curtain Panel 114 cm wide × 70 cm high (optionally as three overlapping panels).
- Double Friction Joint for precise positioning.
- Two integrated Accessory Rails 35 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied); optionally 7.7×25 mm.
- Top-Up Shields included with some products – compatibility see page 16.



312/E-020



312/E-020/3



114 cm wide Product Code	Curtain Panel		Accessory Rail · 35 cm		Top-Up Shields · 57 cm		Top-Up Shields · 27 cm	
	single-piece	three overl.	10×25 mm 9,5×28,6 mm	7.7×25 mm	Lift-off 312/DS/3,32	Fold-down 312/DS/3,33	Lift-off 312/DS/3,30	Lift-off 312/DS/3,41
312/E-020	●		●					
312/E-020/P		●	●					
312/E-020/8	●			●				
312/E-020/8P		●		●				
312/E-020/3	●		●		●		●	●
312/E-020/3/P		●	●		●		●	●
312/E-020/3/8	●			●	●		●	●
312/E-020/3/8P		●		●	●		●	●
312/E-021	●		●		●			
312/E-021/P		●	●		●			
312/E-021/8	●			●	●			
312/E-021/8P		●		●	●			
312/E-022	●		●			●		
312/E-022/P		●	●			●		
312/E-022/8	●			●		●		
312/E-022/8P		●		●		●		

● included · Further top-up shields available – compatibility see page 16

137 cm wide protective shields

Consisting of:


- Curtain Panel 137 cm wide × 75 cm high.
- Double Friction Joint for precise positioning.
- Two integrated Accessory Rails 57 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied).
- Top-Up Shields included with some products – compatibility see page 16.



312/DS-039/4



312/DS-039/4/3



312/DS-039/4

137 cm wide Product Code	Curtain Panel	Accessory Rail · 57 cm	Top-Up Shields · 57 cm		Top-Up Shields · 27 cm	
	single-piece	10×25 mm 9,5×28,6 mm	Lift-off 312/DS/3,32	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30	Lift-off 312/DS/3,41
312/DS-039/4	●	●				
312/DS-039/4/3	●	●		●	●	●
312/DS-0412/4	●	●	●		●	●

● included · Further top-up shields available – compatibility see page 16

162 cm wide protective shields

Consisting of:


- Curtain Panel 162 cm wide × 75 cm high.
- Two Double Friction Joints for precise positioning.
- Two integrated Accessory Rails 57 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied).
- Top-Up Shields included with some products – compatibility see page 16.



312/DS-039/9



312/DS-039/9/4



312/DS-039/9

162 cm wide Product Code	Curtain Panel	Accessory Rail · 57 cm	Top-Up Shields · 57 cm	Top-Up Shields · 27 cm	
	single-piece	10×25 mm 9,5×28,6 mm	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30	Lift-off 312/DS/3,41
312/DS-039/9	●	●			
312/DS-039/9/4	●	●	●	● ×2	●

● included · Further top-up shields available – compatibility see page 16

➤ Table-Mounted Shields Collision-Protected Configurations

Collision-protected lower-body shields have a patented mechanism: The swivel arm deflects upward on collision with the C-arm – so neither C-arm nor shield is damaged and the examination is not interrupted. Each shield is supplied with a flexible 0.5 mm Pb protective curtain.

92 cm · collision-protected

Consisting of:

- Curtain Panel 92 cm wide × 70 cm high (312/DS-008 & /3: 75 cm; optional three overlapping panels).
- Double Friction Joint for precise positioning.
- Patented collision protection – swivel arm deflects upward on collision.
- Top-Up Shields included with some products – compatibility see page 16.



312/E-013



312/DS-008



312/DS-008

92 cm wide Product Code	Curtain Panel		Top-Up Shields · 57 cm		Top-Up Shields · 27 cm
	two-piece	three overl.	Lift-off 312/DS/3,32	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30
312/E-013	●		●		
312/E-013/P		●	●		
312/DS-008	●				
312/DS-008/3	●			●	● ×2

● included · Further top-up shields available – compatibility see page 16

114 cm · collision-protected

Consisting of:


- Curtain Panel 114 cm wide × 70 cm high (optional three overlapping panels).
- Double Friction Joint for precise positioning.
- Two integrated Accessory Rails 35 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied); optional 7.7×25 mm.
- Patented collision protection; Top-Up Shields – compatibility see page 16.



312/E-023



312/E-023/P



312/E-023 (shown without top)

114 cm wide Product Code	Curtain Panel		Accessory Rail · 35 cm		Top-Up Shields · 57 cm
	two-piece	three overl.	10×25 mm 9,5×28,6 mm	7.7×25 mm	Lift-off 312/DS/3,32
312/E-023	●		●		●
312/E-023/P		●	●		●
312/E-023/8	●			●	●
312/E-023/8P		●		●	●


● included · Further top-up shields available – compatibility see page 16

137 cm · collision-protected

Standard version



312/DS-039/6



312/DS-039/6

137 cm wide Product Code	Curtain Panel	Accessory Rail · 57 cm	Top-Up Shields · 57 cm		Top-Up Shields · 27 cm	
	two-piece	10×25 mm 9,5×28,6 mm	Lift-off 312/DS/3,32	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30	Lift-off 312/DS/3,41
312/DS-039/6	●	●				
312/DS-039/6/3	●	●		●	●	●
312/DS-0412/6	●	●	●		●	●

● included · Further top-up shields available – compatibility see page 16

For laterally travelling tables · deep bracket


Deep mounting bracket for transverse table movement. Curtain Panel 137 cm × 75 cm, Double Friction Joint; two Accessory Rails 65 cm.



312/DS-060



lateral table movement



312/DS-060

137 cm · deep bracket Product Code	Curtain Panel	Accessory Rail · 65 cm	Top-Up Shields · 57 cm		Top-Up Shields · 27 cm	
	two-piece	10×25 mm 9,5×28,6 mm	Lift-off 312/DS/3,32	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30	Lift-off 312/DS/3,122
312/DS-060	●	●				
312/DS-060/1/3	●	●		●	●	●

● included · Further top-up shields available – compatibility see page 16

162 cm · collision-protected

Consisting of:


- Curtain Panel 162 cm wide × 75 cm high.
- Two Double Friction Joints for precise positioning.
- Two integrated Accessory Rails 57 cm: 10×25 mm (fitted) / 9.5×28.6 mm (supplied).
- Patented collision protection; Top-Up Shields – compatibility see page 16.



312/DS-039/10



312/DS-039/10/4



312/DS-039/10

162 cm wide Product Code	Curtain Panel	Accessory Rail · 57 cm	Top-Up Shields · 57 cm	Top-Up Shields · 27 cm	
	two-piece	10×25 mm 9,5×28,6 mm	Fold-down/Lift-off 312/DS/3,115	Lift-off 312/DS/3,30	Lift-off 312/DS/3,41
312/DS-039/10	●	●			
312/DS-039/10/4	●	●	●	● ×2	●

● included · Further top-up shields available – compatibility see page 16

Table-Mounted Shields **Top-Up Shields – Options**

KENEX top-up shields are made of the same 0.5 mm Pb protective material as the curtains of the main shield. A flexible plastic material lets them extend 25 cm beyond the table rail and provides additional, higher and wider protection. Their special shape – wider at the top – prevents gaps when pushed together. The compatibility of top-up and main shields can be found on the following pages.

Lift-off top-ups

Solid stainless-steel centring pins hold the top-up shields securely in position. 27 cm wide, shortened top-up shields are made shorter so as not to obstruct control units on the rail.



**312/DS/3.32 ·
Lift-off 57 cm**



**312/DS/3.30 ·
Lift-off 27 cm**



**312/DS/3.41 ·
Lift-off 27 cm
(shortened)**

Fold-down top-ups

57 cm Fold-down top-ups fold away quickly for easier patient access – without removing and repositioning. Fold-down/Lift-off top-ups can additionally be removed and placed on an optional wall bracket (see the Storage section).



Position of the 312/DS-039/6/3 top-ups



312/DS/3.33 · Fold-down 57 cm

Top-Up Sets

312/E-025 · Standard-Height Set

Set of three Lift-off top-ups at standard height.
Height 25 cm above the accessory rail.

Compatible with 312/E-020, /E-021, /E-022, /E-023,
312/DS-039/6 & variants.



312/E-027 · Reduced-Height Set

Set of three Lift-off top-ups at reduced height.
Height 15 cm above the accessory rail.

Compatible with 312/E-020, /E-021, /E-022, /E-023,
312/DS-039/6 & variants.



> Top-Up Shields Heights & Types

Top-up shields are available in standard height (25 cm above the accessory rail) and reduced height (15 cm) – each 57 cm or 27 cm wide. The overview below sorts every top-up shield by height, width and version; the corresponding shield compatibility is shown on the following pages.

Standard height · 25 cm

57 cm wide



312/DS/3.32/3
Lift-off



312/DS/3.32
Lift-off



312/DS/3.33
Fold-down



312/DS/3.115
Fold-down / lift-off

Reduced height · 15 cm

57 cm wide



312/DS/3.106
Lift-off



312/DS/3.128
Lift-off



312/DS/3.139
Fold-down



312/DS/3.167
Fold-down / lift-off

27 cm wide



312/DS/3.30
Lift-off



312/DS/3.41
Lift-off (cropped)



312/DS/3.122
Lift-off (cropped)

27 cm wide



312/DS/3.108
Lift-off



312/DS/3.110
Lift-off (cropped)



312/DS/3.130
Lift-off (cropped)



Top-Up Shields

Compatibility · 57 cm

☒ included in the shield
 ☐ compatible / optional
 ☐ ×2 two mountable

Product code	57 cm wide top-up shields							
	Lift-off (collision)		Lift-off		Fold-down (removable, not lift-off)		Fold-down / lift-off	
	Standard height 312/DS/3.32/3	Reduced height 312/DS/3.106	Standard height 312/DS/3.32	Reduced height 312/DS/3.128	Standard height 312/DS/3.33	Reduced height 312/DS/3.139	Standard height 312/DS/3.115	Reduced height 312/DS/3.167
78 cm wide								
312/E-005			☐	☐	☐	☐		
312/E-006			☒	☐	☐	☐		
312/E-007			☐	☐	☒	☐		
92 cm wide								
312/E-010			☐	☐	☐	☐		
312/E-011			☒	☐	☐	☐		
312/E-012			☐	☐	☒	☐		
100 cm wide								
312/E-015			☐	☐	☐	☐		
312/E-016			☒	☐	☐	☐		
312/E-017			☐	☐	☒	☐		
114 cm wide								
312/E-020			☐	☐	☐	☐		
312/E-021			☒	☐	☐	☐		
312/E-022			☐	☐	☒	☐		
137 cm wide								
312/DS-039/4			☐	☐			☐	☐
312/DS-039/4/3			☐	☐			☒	☐
312/DS-0412/4			☒	☐			☐	☐
162 cm wide								
312/DS-039/9			☐	☐			☐	☐
312/DS-039/9/4			☐	☐			☒	☐
92 cm wide · collision-protected								
312/E-013			☒	☐			☐	☐
312/DS-008			☐	☐			☐	☐
312/DS-008/3			☐	☐			☒	☐
114 cm wide · collision-protected								
312/E-023			☒	☐			☐	☐
137 cm wide · collision-protected								
312/DS-039/6			☐	☐			☐	☐
312/DS-039/6/3			☐	☐			☒	☐
312/DS-0412/6			☒	☐			☐	☐
312/DS-060			☐	☐			☐	☐
312/DS-060/1/3			☐	☐			☒	☐
162 cm wide · collision-protected								
312/DS-039/10			☐	☐			☐	☐
312/DS-039/10/4			☐	☐			☒	☐



Top-Up Shields

Compatibility · 27 cm

☒ included in the shield
 ☐ compatible / optional
 ☐ ×2 two mountable

Product code	27 cm wide top-up shields					
	Lift-off		Lift-off (cropped)		Lift-off (cropped)	
	Standard height 312/DS/3.30	Reduced height 312/DS/3.108	Standard height 312/DS/3.41	Reduced height 312/DS/3.110	Standard height 312/DS/3.122	Reduced height 312/DS/3.130
78 cm wide						
312/E-005	☐	☐				
312/E-006	☐	☐				
312/E-007	☐	☐				
92 cm wide						
312/E-010	☐ ×2	☐ ×2				
312/E-011	☐ ×2	☐ ×2				
312/E-012	☐ ×2	☐ ×2				
100 cm wide						
312/E-015			☐	☐		
312/E-016			☐	☐		
312/E-017			☐	☐		
114 cm wide						
312/E-020	☐	☐	☐	☐		
312/E-021	☐	☐	☐	☐		
312/E-022	☐	☐	☐	☐		
137 cm wide						
312/DS-039/4	☐	☐	☐	☐		
312/DS-039/4/3	☒	☐	☒	☐		
312/DS-0412/4	☒	☐	☒	☐		
162 cm wide						
312/DS-039/9	☐ ×2	☐ ×2	☐	☐		
312/DS-039/9/4	☒ ×2	☐ ×2	☒	☐		
92 cm wide · collision-protected						
312/E-013	☐ ×2	☐ ×2				
312/DS-008	☐ ×2	☐ ×2				
312/DS-008/3	☒ ×2	☐ ×2				
114 cm wide · collision-protected						
312/E-023	☐	☐	☐	☐		
137 cm wide · collision-protected						
312/DS-039/6	☐	☐	☐	☐		
312/DS-039/6/3	☒	☐	☒	☐		
312/DS-0412/6	☒	☐	☒	☐		
312/DS-060	☐	☐			☐	☐
312/DS-060/1/3	☒	☐			☒	☐
162 cm wide · collision-protected						
312/DS-039/10	☐ ×2	☐ ×2	☐	☐		
312/DS-039/10/4	☒ ×2	☐ ×2	☒	☐		



Head-End Shields

Head-end shields significantly reduce the scattered radiation emanating from the table or the patient's body. They consist of a lightweight, durable structure with a flexible protective curtain with 0.5 mm Pb lead equivalent.

- The shield slides onto the carbon table top and holds its position even with a tilted table.
- The swivelling side wings can be shaped around the user, extended forward individually or folded away to save space; the swivel arms can be quickly repositioned for full C-arm movement.
- Each shield is made for a specific table and is optionally available with wall storage.

313/A2 · Custom-made head-end shields

Individual solutions depending on size and table type. Note: A measurement sheet must be completed when ordering (available as a PDF on the product page or on request).

Product code	Wall Bracket
313/A2 · measurement sheet required	☐

☒ included · ☐ optional

313/A2-001 & 313/A2-008 · for Siemens Healthineers Artis zee

313/A2-001 · 120 cm wide

1× 56 cm fixed centre section · 2× 32 cm swivel arms · 80 cm curtain · wall bracket



313/A2-001



folded · wall bracket

313/A2-008 · 122 cm wide

1× 58 cm fixed centre section · 2× 32 cm swivel arms · 80 cm curtain · wall bracket



Product code	Wall Bracket
313/A2-001 · for 525 mm Artis / 520 mm TruSystem 7500 (Artis zee)	☒
313/A2-008 · for 550 mm Artis zee Multipurpose (MP)	☒

☒ included · ☐ optional

313/A2-010 · for GE Healthcare Elegance


Consisting of: 1× table adapter bracket · 1× 37 cm fixed centre section · 2× swivel arms (24 & 68 cm wide) · 60 cm curtain · 3× 29 cm Lift-off top-ups · wall bracket.



313/A2-010 · without top-ups



313/A2-010 · wall bracket



Product code	Wall Bracket	29 cm Lift-off top-ups
313/A2-010 · for GE Healthcare Elegance	●	● ×3

● included · ○ optional

313/A2-014 · for Canon Medical (Toshiba) CAT-880B

114 cm wide · 1× 50 cm fixed centre section · 2× 32 cm swivel arms · 78.5 cm curtain · wall bracket.



Product code	Wall Bracket
313/A2-014 · for 450 mm tilt/cradle table CAT-880B	☒

☒ included · ☐ optional

313/A2-018 · for Baxter Carbon FloatLine & SQ14 X-TRA

120 cm wide · 1× 56 cm fixed centre section · 2× 32 cm swivel arms · 80 cm curtain · wall bracket.



Product code	Wall Bracket
313/A2-018 · for 520 mm Carbon FloatLine & SQ14 X-TRA	☒

☒ included · ☐ optional

313/A3 · for Philips Allura AD7

Consisting of: 121 cm wide · 1× 57 cm fixed centre section · 2× 32 cm swivel arms · 80 cm curtain · 1× 57 cm Fold-down top-up.



Product code	Wall Bracket	57 cm Fold-down top-up
313/A3-001 · for 500 mm Allura AD7	☐	☒
313/A3-002 · for 500 mm Allura AD7 (incl. wall bracket)	☒	☒

☒ included · ☐ optional



Foot-End Shield

The upper foot-end shield protects against scattered radiation and additionally against possible blood splashes during interventional procedures.

- For use across the table to the foot end.
- Fits both 25 mm high table rails (7.7–10 mm thick).
- Continuous lead equivalent: 0.5 mm Pb.

Note: A measurement sheet must be completed when ordering (available as a PDF on the product page or on request).

311/TC/004 · Upper foot-end table shield

1× Lead-acrylic window approx. 60 cm high × 80 cm wide · 1× flexible protective curtain below the acrylic, 20 cm high · 2× handles for positioning.



311/TC/004



311/TC/004

> Table-Mounted Shields Storage · Wall Bracket

312/W-LBT · Wall-bracket set

Wall-bracket set, consisting of:

- 1× 312/W-LB (for lower-body shield)
- 1× 312/W-TOP (for top-up shields)




312/W-LBT

312/W-LB · Wall bracket for lower-body shield

Practical storage location when shields are removed from the table. Holder for 1× KENEX lower-body shield (type 312).




312/W-LB



312/W-TOP · Wall bracket for top-up shields

Storage location for removed top-up shields. Holder for 1× 57 cm + 2× 27 cm top-up shields.



Compatibility · 57 cm top-up shields

Product Code	312/DS/3.32/3	312/DS/3.106	312/DS/3.32	312/DS/3.128	312/DS/3.33	312/DS/3.139	312/DS/3.115	312/DS/3.167
312/W-TOP	○	○	○	○			○	○

Compatibility · 27 cm top-up shields

Product Code	312/DS/3.30	312/DS/3.108	312/DS/3.41	312/DS/3.110	312/DS/3.122	312/DS/3.130
312/W-TOP	○	○	○	○	○	○

○ compatible



Storage · Mobile Stands

Mobile stands for departments with little or unsuitable space for wall brackets. Small footprint when parked, large lockable castors with front guard against floor cables as well as a heavy, tip-resistant base.

1342/A2/E-001

1× table shield (Type 312) · 1× top-up set (1×57 + 2×27 cm) · 1× head-end shield (313/A2-008 or 313/A2-001)



1342/A2/E-001

1342/A3/E-001

1× table shield (Type 312) · 1× top-up set (1×57 + 2×27 cm) · 1× head-end shield (313/A3-001)



1342/A3/E-001

1342/A2/F

2× lower-body shields (Type 312) · 2× top-up sets (2×57 + 4×27 cm)



1342/A2/F



Product Code	Lower-body shield Type 312	Head-end shield Type 313	Top-up set 1×57 + 2×27 cm
1342/A2/E-001 · compatible with head-end shield 313/A2-008 & 313/A2-001	○	○	○
1342/A3/E-001 · compatible with head-end shield 313/A3-001	○	○	○
1342/A2/F	○ ×2		○ ×2 sets

○ compatible

CONTACT

Your Partner for Radiation Protection.

We source and supply certified medical devices. Beyond that, we procure almost any product independently of the manufacturer and equip entire areas on a project basis.

medeqon GmbH

Engineering Office for Medical Technology

Bergstrasse 42/5/3

2102 Hagenbrunn · Austria

office@medeqon.com

+43 1 3580045

www.medeqon.com